


Listing of Equivalent Items available from **CHOICE ANALYTICAL** for use on VARIAN instrumentation



Sundry Items used in Atomic Spectroscopy Instruments ICP, AAS etc. and associated Accessories and Auto Samplers etc.

VARIAN Part No.	Part No. CHOICE	Description	Unit of Sale
-----------------	-----------------	-------------	--------------

Sampling and General Use Tubing

—	4061-01	Teflon Tubing, 3/16 id x 5/16 od	per metre
—	4097-302	Teflon Tubing, .026 id x .050 od	each
—	4097-322	Teflon Tubing, .063 id x .125 od	per metre
—	4060-38TT	Teflon Tubing, .6mm id, 10ft	each

Deuterium Lamp & Photomultiplier Lamp Arsenic

56-101215-00	P-706	D2 Lamp Spectra AA & x-75	each
61-901761-00	4085-08	Photomultiplier Tube, Arsenic	each
20-100150-00	4097-16	Quartz Window	each

Sample Cups / Tubes used on Auto Samplers

VARIAN Part No.	Part No. CHOICE	Description	Unit of Sale
	AAA-1001	2ml Polystyrene Conical Cups	pkt/1,000
	AAA-1010	2ml Polystyrene Conical Cups	pkt/10,000
99-100282-00	AAA-2001	2ml Polyethylene Conical Cups	pkt/1,000
99-100282-00	AAA-2010	2ml Polyethylene Conical Cups	pkt/10,000
99-100284-00	AAA-4001	4ml Polystyrene Conical Cups	pkt/1,000
—	4092-320	16-17mm x 100mm, Polypropylene Tube	pkt/1,000
—	4092-321	16-17mm x 100mm, Polypropylene Tube	pkt/10,000
—	4092-324	13 x 100mm Polystyrene Tube	pkt/1,000
—	4092-327	13 x 100mm Polypropylene Tube	pkt/1,000
—	4092-330	12 x 75mm Polypropylene Tube	pkt/1,000
—	4092-400	Centrifuge Tube, 50ml, Polypropylene, metal free	pkt/500
—	SC0045	4.5ml Polystyrene Tube (10 x 65mm)	pkt/500
—	SC0060	13 x 100mm Polypropylene Tube 6ml	pkt/1,000
		16 x 125mm Polypropylene Tube 16ml	pkt/1,000

Note: For higher acidic or organic samples **Polyethylene** is recommended whereas for low acidity aqueous samples **Polystyrene** is suggested.

Peristaltic Pump Tubing:-

- Tygon / PVC
- Coated PVC
- Viton
- Silicone

Suitable for use on:-

VARIAN

**AAS and VGA,
ICP-OES, ICP-MS Instruments**

All tubing features **two or three colour-coded stops** for easy determination of size and is available in • **Standard Tygon / PVC** • **Solvent Resistant coated PVC (yellow tinted)** • **Viton (black opaque)** • **Silicone**. Tubing overall length is 16", with a 6" centre pumping section separated by two colour-coded stops, three stops (72 mm between stops) Our Viton Tubing (black opaque) has an overall length of 7" and also includes two stops. These products are guaranteed for quality and normal life.

Tab/Stop Colour Code	ID inches/mm	Part No. VARIAN	Standard Tygon/PVC Tubing	Coated PVC Solvent Resistant Tubing	Viton High Chem. & Temperature Resistant	Pack Size	Silicone Tubing	Pack Size
TWO-STOP-TYPE								
Black/Black	0.030/0.76	37-100272-00	AR-4308	–	–	12/pk	–	–
Black/Black	0.030/0.76	99-100610-00	–	–	AR-4308B	12/pk	–	–
Black/Black	0.030/0.76	37-100348-00	–	AR-4308Y	–	12/pk	AR-4308S	6/pk
White/White	0.040/1.02	37-100344-00	AR-4311	–	–	12/pk	–	–
White/White	0.040/1.02	37-100350-00	–	–	AR-4311B	12/pk	–	–
White/White	0.040/1.02	37-100350-00	–	AR-4311Y	–	12/pk	AR-4311S	6/pk
Grey/Grey	0.051/1.30	37-100345-00	AR-4315	–	–	12/pk	–	–
Grey/Grey	0.051/1.30	37-100352-00	–	–	AR-4315B	12/pk	–	–
Grey/Grey	0.051/1.30	37-100352-00	–	AR-4315Y	–	12/pk	AR-4315S	6/pk
Blue/Blue	0.065/1.65	37-100346-00	AR-4318	–	–	12/pk	–	–
Blue/Blue	0.065/1.65	99-100609-00	–	–	AR-4318B	12/pk	–	–
Purple/Purple	0.081/2.06	–	AR-4321	–	–	12/pk	–	–
Purple/Black	0.090/2.29	37-100273-00	AR-4322	–	–	12/pk	–	–

Many other sizes and types of Tubing and Pump Tubes are also available e.g. **3 stop type**. Please enquire providing full details of requirement.

