

Listing of Equivalent Items  
available from **CHOICE ANALYTICAL** for use on:-



ELTRA CS 800

## Consumables Listing for ELTRA CS 800 Carbon / Sulphur Analyser

In sequence of

ELTRA Part No.	CHOICE Part No.	Description	Unit of Sale
<b><u>FURNACE CLEANING ASSEMBLY</u></b>			
14021	EL-14021	Upper Furnace Block	each
14045	EL-14045	Brush	each
14051	EL-14051	Brush Holder	each
14072	EL-14072	Ceramic Heat Shield for Brush	each
14080	EL-14080	Cleaning Mechanism Rod	each
70120	EL-70120	O Ring	each
71031	EL-71031	Wire Brush for cleaning 14072 Heat Protector	each
75120	EL-75120	Spring	each
75130	EL-75130	Safety Spring	each
75150	EL-75150	Metal Tube	each

### **FURNACE**

11105	EL-11105	Metal Dust Filter	each
13067	EL-13067	Combustion Coil	each
14009	EL-14009	Pneumatic Cylinder for Furnace Lift	each
14021	EL-14021	Upper Furnace Block	each
14025	EL-14025	Lower Furnace Block	each
14090	EL-14090	Bearing	each
14100	EL-14100	Mounting Rod	each
14130	EL-14130	Combustion Tube	each
14168	EL-14168	Pedestal	each
14170	EL-14170	Pedestal Mount	each
14180	EL-14180	Furnace Closure	each
14185	EL-14185	Tray	each
14200	EL-14200	Metal Tube	each
14210	EL-14210	Threaded Rod	each
70380	EL-70380	O Ring	each
70390	EL-70390	O Ring	each
77210	EL-77210	Oscillator Tube	each
90150	AR-3818	Ceramic Crucibles (1000)	1000/Box
	AR-040	Ceramic Lids 4mm hole (250)	250/Box
	AR-043	Ceramic Lids 10mm hole (250)	250/Box

... Ctd. Page 2

## Consumables Listing for ELTRA CS 800 Carbon / Sulphur Analyser

In sequence of ELTRA Part No.	CHOICE Part No.	Description	Unit of Sale
<b>COMBUSTION ACCELERATORS &amp; CHEMICALS</b>			
90200	CA-171	Anhydrone (454gm)	454gm
	CA-008	Tungsten II Accelerator (2.5kg)	2.5kg
90210	CA-2174	CO2 Absorber (500gm) CA-sorb A (14-25 mesh)	500gm
90220	CA-027	Tungsten Accelerator (2.5kg)	2.5kg
90260	CA-077	Iron Chip Accelerator (0.9kg)	0.9kg
	AR-098	Cellouse Wool (454gm)	454gm
	AR-170	Rare Earth Copper Oxide (50gm)	50gm
90280	AR-076	Tin Accelerator (908gm - 2lb)	908gm
90330	AR-177-50	Quartz Wool (50gm)	50gm
	AR-263	Copper Accelerator (1.35kg)	1.35kg
	AR-579	Metal Accelerator Scoop	each
	AR-929	Crucible Tongs	each
90331	CA-081	Glass Wool (454gm)	454gm

<b>DETERMINATOR</b>			
11064	EL-11064	Reagent Tube	each
11091	EL-11091	Dust Trap Mechanism	each
11105	EL-11105	Metal Dust Trap	each
11110	EL-11110	Filter Housing	each
11115	EL-11115	Furnace Outlet Tube	each
11120	EL-11120	Paper Filter Holder	each
11185	EL-11185	Paper Filter (10)	10/Pack
11480	EL-11480	Adjustable Restrictor	each
15083	EL-15083	Gas Flow Indicator 151/h	each
15087	EL-15087	Gas Flow Indicator 300/h	each
20000	EL-20000	Catalyst Furnace	each
20040	EL-20040	Catalyst Tube	each
70210	EL-70210	O Ring	each
70230	EL-70230	O Ring	each
70320	EL-70320	O Ring	each
70370	EL-70370	O Ring	each
71030	EL-71030	Soft Brush for Cleaning 11105 Metal Dust Trap	each
72010	EL-72010	Pressure Gauge	each
72020	EL-72020	Pressure Gauge	each
77430	EL-77430	Panel Meter 100uA	each
78010	EL-78010	Mains Switch	each

