


Listing of Equivalent Items  
available from **CHOICE ANALYTICAL** for use on:-



LECO S/SC-144

### Consumables Listing for LECO S-144, SC-144 and S-144DR

LECO Part No.	Part No. CHOICE	Description	Unit of Sale
	ASCRM 012A-2	Coal Standard 0.33%S	125gm
	ASCRM 012B-2	Coal Standard 0.66%S	125gm
	ASCRM 012C-2	Coal Standard 1.11%S	125gm
	ASCRM 012D-2	Coal Standard 5.21%S	125gm
	ASCRM 013-9	Coal Standard 0.59%S + other results	250gm
175-590	AR-5590	Pump Repair Kit	each
501-001	AR-1701	Coal Standard 0.4%S	50gm
501-005	AR-1704	Coal Standard 1.0%S	50gm
501-020	AR-1708	Coal Standard 3.0%S	50gm
501-034-500	AR-034-500	Standard 12.0%C	500gm
501-062	AR-061	Boat Puller (Manual)	each
501-078	AR-078	Iron Powder Accelerator	454gm = 1lb
501-081	CA-081	Glass Wool	454gm = 1lb
501-171	CA-171	Anhydrone 8-20# Mesh	454gm
501-241	AR-241	Silicone Grease	5.3 oz
501-426	CA-426	Comb. Control (coarse) for Solids	454gm
501-427	CA-427	Comb. Control (fine) for Liquids	454gm
501-614	AEE-2642	Spatula	each
501-636	AR-636	Vanadium Pentoxide	250gm
502-023	AR-473	Reagent Funnel	each
502-085	AR-3085	Zinc Sulphide Std.~30%S	50gm
502-301	AR-2301	Tungsten Oxide	125gm
502-321	AR-321	Alphacat Accelerator	250gm
502-343	CA-2343	Nickel Boat Liners for 528-203 (AR-8203) Boat	100/Pk
503-032	AR-901	Glass Scoop	each
529-204	AR-9204P	Combustion Boat	(Bag of 50)
529-204-250	AR-9204	Combustion Boat (note AR-9204 is Box of 500)	(Box of 500)
537-103	AR-103	Polyflo Tubing 1/4" OD	per foot
601-483	AR-483	Printer Paper, Sheets	600
601-484	AR-484	Printer Ribbon	each
602-615	AR-2615	Ceramic Furnace Housing	each
606-224	AR-6224	Thermocouple	each
606-286	AR-6286	Short Strap	each
606-319	AR-6319	Lead Wire Braided 6.5	each
606-320	AR-6320	Clamp (Furn. Element)	each
606-327	AR-6327	Gas Dessicator Tube	each

ctd. page 2



**CHOICE Analytical Pty Ltd**

ABN 66 052 203 719

Supplies for Laboratory Instruments

Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : [info@choiceanalytical.com.au](mailto:info@choiceanalytical.com.au)Web : [www.choiceanalytical.com.au](http://www.choiceanalytical.com.au)

## Consumables Listing for LECO S-144, SC-144 and S-144DR

LECO Part No.	Part No. CHOICE	Description	Unit of Sale
606-333	AR-6333	O-Ring	each
606-340	AR-6340	Lance, 250mm length, 1 hole	each
606-403	AR-403	Brush	each
606-601	AR-6601	Heating Element (set)	set of 4
606-602	AR-6602	End Cap for Ceramic Furnace Housing	each
606-655	AR-6655	Large Gas Drying Tube	each
611-476	AR-1476	O-Ring	each
611-477	AR-1477	O-Ring	each
616-068	AR-6068	O-Ring, Inner Combustion Tube S/SC-144	each
616-070	AR-6070	Inner Combustion Tube S/SC-144	each
616-086	AR-6086	Stepped Boat Stop, S/SC-144	each
616-089	AR-6089	Outer Combustion Tube, S/SC-144	each
760-193	AR-193	Screen for reagent tube	each
761-747	AR-611	Schutze's Reagent	100gm
761-929	AR-929	Tongs Boat pickup	each
763-265 see 501-081	CA-081	Glass wool (763-265=1/2lb) Sold only as CA-081 (501-08=1lb)	454gm=1lb
768-980	AR-980	Particle Filter	each
769-608	AR-608	Antimony Halogen Trap Material	200gm
769-610	AR-610	Halogen Trap Material	225gm
772-520	AR-520	O-Ring	each
773-320	AR-320	O-Ring	each
776-978	AR-978	Large Glass Scoop	each
781-273	AR-028	Alpha Plate (for hot items)	each
781-335	AR-1335	Quartz Combustion Boat	each
782-058	now 502-343	Nickel Liners for 528-203 (AR-8203) Boat	100/pk
783-897	AR-3897	O-Ring Outer Combustion Tube	each

### Heating Elements



**AR-6318** (606-318) Set of 3 Elements

**AR-6601** (606-601) Set of 4 Elements


**AR-6320** (606-320) Element Braid Clamp (each)

**AR-6286** (606-286) Short Strap (each)

### AR-2615 (602-615) Ceramic Furnace Housing.

For use on the SC432 and SC444.

This Furnace Housing is manufactured from high purity ceramic material rated at 1600°C.



### AR-6602 (606-602) End Cap for Ceramic Furnace Housing.

Made from the same high purity ceramic material as the Housing.

AR-6602

AR-6602

AR-2615



**CHOICE Analytical Pty Ltd**

Supplies for Laboratory Instruments

Tel : + 61 2 9980 6240

Fax : + 61 2 9980 6344

E-mail : [info@choiceanalytical.com.au](mailto:info@choiceanalytical.com.au)

Web : [www.choiceanalytical.com.au](http://www.choiceanalytical.com.au)

ABN 66 052 203 719