



LOW FORM CRUCIBLES (TRANSLUCENT or * TRANSPARENT)

Part No.*	Description	Type	EXTERNAL DIMENSIONS		
			Height (H)	Diameter (D)	Capacity
S 1410	Low form crucibles	MI/LFC/05	19 mm	29 mm	5 ml
S 1420	Low form crucibles	MI/LFC/10	19 mm	38 mm	10 ml
S 1430	Low form crucibles	MI/LFC/15	25 mm	41 mm	15 ml
S 1440	Low form crucibles	MI/LFC/25	28 mm	47 mm	25 ml
S 1450	Low form crucibles	MI/LFC/50	37 mm	57 mm	50 ml
S 1460	Low form crucibles	MI/LFC/80	45 mm	67 mm	80 ml
S 1470	Low form crucibles	MI/LFC/150	51 mm	82 mm	150 ml

*If required in TRANSPARENT use 'T' instead of 'S'.



LIDS FOR CRUCIBLES

Part No.*	Description	Type
S 1415	Lid for 5 ml crucible	MI/LFC/05-L
S 1425	Lid for 10 ml crucible	MI/LFC/10-L
S 1435	Lid for 15 ml crucible	MI/LFC/15-L
S 1445	Lid for 25 ml crucible	MI/LFC/25-L
S 1455	Lid for 50 ml crucible	MI/LFC/50-L
S 1465	Lid for 80 ml crucible	MI/LFC/80-L
S 1475	Lid for 150 ml crucible	MI/LFC/150-L

*If required in TRANSPARENT use 'T' instead of 'S'.

Other configurations, shapes, sizes, translucent, transparent, glazed, unglazed etc. are available – please enquire and detail your specific requirements.

