


90-5655-xx
Platinum Crucible



90-5665-xx
Platinum Cover



90-5675-xx 90-5677xx
Reshapers

For Determination of Volatile Matter in COAL ASTM D-3175-89a (Reapproved 1993)

Platinum Crucibles and Covers

Part No.	Description	Capacity (mls)	Diameter (mm)	Height (mm)	Approx Weight (gms)
90-5655-34	Pt Crucible	15	30	31	11
90-5665-34	Pt Cover (for above)				5
90-5655-39	Pt Crucible	20	33	35	15
90-5665-39	Pt Cover (for above)				7

Above items conform to ASTM D-3175 6.Apparatus.

6.1 Platinum Crucible, with closely fitting cover, for coal. The crucible shall be of not less than 10 or more than 20-ml capacity, not less than 25 or more than 35 mm in diameter, and not less than 30 or more than 35 mm in height.

Note: Other platinum crucibles and covers (10, 25 and 30 ml capacity) are also available if required but fall outside above specification.

RESHAPERS - For Platinum Crucibles and Insert Covers

Durable, plastic reshapers are available for all sizes of crucibles and insert covers listed under 90-5655-xx and 90-5665-xx. They remove dents and distortions by means of simple pressure in the mold. Reshapers are designed for an exact fit with the vessel to be reshaped. As a result, it is best to order the reshaper when ordering crucibles and covers. When ordering reshapers needed for your existing platinum ware, specify the complete catalog number of your crucibles or covers.

RESHAPER MOLDS - for above Platinum Crucibles and Covers

Reshaper Part No.	For use with Crucible Type	For use with Cover Type	Unit of Sale
90-5675-34	90-5655-34		each
90-5675-39	90-5655-39		each
90-5677-34		90-5665-34	each
90-5677-39		90-5665-39	each