


### Peristaltic Pump Tubing:-

- Tygon/PVC
- Coated PVC
- Viton
- Silicone

Suitable for use on:-

**VARIAN**

**AAS and VGA,  
ICP-OES, ICP-MS Instruments**



All tubing features **two or three colour-coded stops** for easy determination of size and is available in • **Standard Tygon/PVC** • **Solvent Resistant coated PVC (yellow tinted)** • **Viton (black opaque)** • **Santoprene (white opaque)** • **Silicone**. Tubing overall length is 16", with a 6" centre pumping section separated by two colour-coded stops. Our Viton Tubing (black opaque) has an overall length of 7" and also includes two stops. These products are guaranteed for quality and normal life.

Tab/Stop Colour Code	ID Inches/mm	Varian Part No.	Standard Tygon/PVC Pack of 12	Coated PVC Solvent Tubing Pack of 12	Viton High Chemical Resistant Pack of 12	Santoprene Tubing Pack of 12	Silicone Tubing Pack of 6
Orange/Yellow	0.020/0.51	37-100435-00	-	-	116-0651-05	-	-
Orange/Yellow	0.020/0.51	37-100440-00	-	-	-	SC0322	-
Black/Black	0.030/0.76	37-100272-00	116-0549-07	-	-	-	-
Black/Black	0.030/0.76	37-100436-00	-	-	116-0651-07	-	-
Black/Black	0.030/0.76	37-100348-00	-	116-0533-07	-	-	116-0497-08
Black/Black	0.030/0.76	37-100441-00	-	-	-	SC0324	-
White/White	0.040/1.02	37-100344-00	116-0549-09	-	-	-	-
White/White	0.040/1.02	37-100437-00	-	-	116-0651-09	-	-
White/White	0.040/1.02	37-100350-00	-	116-0533-09	-	-	116-0497-10
White/White	0.040/1.02	37-100442-00	-	-	-	SC0326	-
White/White	0.040/1.02	99-100609-00	-	-	-	SC0326	-
Grey/Grey	0.051/1.30	37-100345-00	116-0549-11	-	-	-	-
Grey/Grey	0.051/1.30	37-100438-00	-	-	116-0651-11	-	-
Grey/Grey	0.051/1.30	37-100352-00	-	116-0533-11	-	-	116-0497-12
Grey/Grey	0.051/1.30	37-100443-00	-	-	-	SC0328	-
Grey/Grey	0.051/1.30	99-100610-00	-	-	-	SC0328	-
Blue/Blue	0.065/1.65	37-100346-00	116-0549-13	-	-	-	-
Blue/Blue	0.065/1.65	37-100439-00	-	-	116-0651-13	-	-
Blue/Blue	0.065/1.65	37-100444-00	-	-	-	SC0331	-
Purple/Purple	0.081/2.06	-	116-0549-15	-	-	-	-
Purple/Black	0.090/2.29	37-100273-00	116-0549-16	-	-	-	-

Many other sizes of Tubing and Pump Tubes are available.