


Peristaltic Pump Tubing:-

- Tygon / PVC
- Coated PVC
- Viton
- Silicone

Suitable for use on:-

- Auto-Analyzers
- AAS & ICP Instruments etc.



All tubing features **three colour-coded stops** for easy determination of size and is available in • **Standard Tygon / PVC** • **Solvent Resistant coated PVC (yellow tinted)** • **Viton (black opaque)** • **Silicone** • **Santoprene [*] (creamy opaque)**. Tubing overall length is 18", with two 3" centre pumping sections separated by three colour-coded stops. Our Viton Tubing (black opaque) has an overall length of 7" and also includes three stops. These products are guaranteed for quality and normal life.

Tab/Stop Colour Code	ID Inches/mm	Standard Tygon/PVC Pack of 12	Coated PVC Solvent Resistant Tubing Pack of 12	Viton High Chemical & Temperature Resistant Pack of 12	Santoprene Tubing * Pack of 12	Silicone Tubing Pack of 6
Orange/Black	0.005/0.13	AR-34300	AR-34300Y	-	-	-
Orange/Red	0.008/0.19	AR-34301	AR-34301Y	-	-	-
Orange/Blue	0.010/0.25	AR-34302	AR-34302Y	-	AR-34302SA/12	-
Orange/Green	0.015/0.38	AR-34303	AR-34303Y	-	AR-34303SA/12	-
Green/Yellow	0.017/0.44	AR-34304	AR-34304Y	-	AR-34304SA/12	--
Orange/Yellow	0.020/0.51	AR-34305	AR-34305Y	AR-34305B	AR-34305SA/12	-
White/Yellow	0.022/0.57	AR-34306	AR-34306Y	-	-	-
Orange/White	0.025/0.64	AR-34307	AR-34307Y	AR-34307B	AR-34307SA/12	AR-34307S
Black/Black	0.030/0.76	AR-34308	AR-34308Y	AR-34308B	AR-34308SA/12	AR-34308S
Orange/Orange	0.035/0.89	AR-34309	AR-34309Y	AR-34309B	AR-34309SA/12	AR-34309S
White/Black	0.037/0.95	AR-34310	AR-34310Y	-	-	-
White/White	0.040/1.02	AR-34311	AR-34311Y	AR-34311B	AR-34311SA/12	AR-34311S
White/Red	0.043/1.09	AR-34312	AR-34312Y	-	-	-
Red/Red	0.045/1.14	AR-34313	AR-34313Y	AR-34313B	AR-34313SA/12	AR-34313S
Red/Grey	0.048/1.22	AR-34314	AR-34314Y	-	-	-
Grey/Grey	0.051/1.30	AR-34315	AR-34315Y	AR-34315B	AR-34315SA/12	AR-34315S
Yellow/Yellow	0.056/1.42	AR-34316	AR-34316Y	AR-34316B	AR-34316SA/12	AR-34316S
Yellow/Blue	0.060/1.52	AR-34317	AR-34317Y	AR-34317B	AR-34317SA/12	AR-34317S
Blue/Blue	0.065/1.65	AR-34318	AR-34318Y	AR-34318B	AR-34318SA/12	AR-34318S
Blue/Green	0.069/1.75	AR-34319	AR-34319Y	-	-	-
Green/Green	0.073/1.85	AR-34320	AR-34320Y	AR-34320B	AR-34320SA/12	AR-34320S
Purple/Purple	0.081/2.06	AR-34321	AR-34321Y	AR-34321B	AR-34321SA/12	AR-34321S
Purple/Black	0.090/2.29	AR-34322	AR-34322Y	AR-34322B	AR-34322SA/12	AR-34322S
Purple/Orange	0.100/2.54	AR-34323	AR-34323Y	AR-34323B	AR-34323SA/12	AR-34323S
Purple/White	0.110/2.79	AR-34324	AR-34324Y	AR-34324B	AR-34324SA/12	AR-34324S
Black/White	0.125/3.17	AR-34325	AR-34325Y	-	-	-

Many other sizes and types of Tubing and Pump Tubes are also available. Please enquire providing full details of requirement

* Santoprene = PharMed = Ismaprene (All are same material - creamy opaque)

